

# AdvanceSeal™ 2000

Effective: August 15, 2025  
Document No.: MFG-PS-457-A

## Dimensional Information and Tolerances

CRITERIA	DIMENSIONAL INFORMATION / SPECIFICATIONS
Roll Stock:	
Portrait Roll (W X L)	77.5mm x 610m ± 0.5 (3.05" x 2,001')
Landscape Roll (W x L)	114.5mm x 500m ± 0.5 (4.51" x 1640')
Core Inner Diameter	3.038"

## Construction

LAYER	DESCRIPTION	THICKNESS
Top coat	N/A	N/A
Film	Coated Aluminum Foil	1mil ±15%
Sealant layer	Polyolefin	1.5mil ±15%

## Activation Temperature Table

PLATE MATERIAL	RECOMMENDED TEMPERATURE	RECOMMENDED TIME
Polypropylene	145°C – 160°C	2 - 3 Seconds
Dual Component (polycarbonate/polypropylene)	145°C – 160°C	2 - 3 Seconds
COC	145°C – 160°C	2 - 3 Seconds
Polystyrene	75°C – 120°C	2 - 3 Seconds

## Applications

- PCR
- Sample storage and transport
- Long term storage
- High throughput applications
- Automated applications
- Anaerobic Conditions
- Re-sealing pierced plates

## Features & Benefits

- Excellent protection from light, gas, and moisture
- Supports resealing with additional layers of identical foil
- Strong sample protection, even for small sample volumes
- Pierceable aluminum foil
- Recommended for use at temperatures from -80°C to +105°C
- Certified free of detectable levels of RNase, DNase, and Human DNA (PCR Certified)
- PCR inhibitor free
- Endotoxin free
- Excellent chemical resistance