

1x8 MiniStrips™

For sealing only one 8-well row of a 96-well plate



Application:

General laboratory use, incubation and ELISA



Temperature Range:

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
GP-CS1X8-400	1x8 MiniStrips	400 Strips
GPS-CS1X8-400	1x8 MiniStrips, Sterile	400 Strips



Features and Benefits

- Certified DNase- & RNase- free
- A pliable polypropylene film that minimizes aerosols upon removal
- Each sheet contains 8 strips on a continuous liner for ease of use
- Apply individually or two-at-a-time
- Less expensive than full-size sheets when using only a partial plate
- Designed for use with either 1x8 strip-well ELISA plates or standard 96-well plates where rows must be selectively protected or accessed
- Use as a secondary film to seal previously pierced wells or to prevent accidental pipetting
- Ideal for COVID-19 testing to seal a row of 8 wells after swabs are deposited into a deep-well plate

Bundling Opportunities

- 96-well plates
- ELISA plates
- ELISA strip-well plates
- Plate washers
- Plate readers
- Pipettors and pipet tips
- Reservoirs