

Multi-Swipe™ Sealing Paddle

For sealing flat-top and raised-rim plates

i Application:
Microplate sealing

Material:
Polyethylene

CATALOG #	DESCRIPTION	QUANTITY
PDL-5PK	Paddle, 1 bag of 5	5 Paddles



Features and Benefits

- Reusable, recyclable
- Ensures a secure and uniform, pressure-sensitive seal around all wells
- Designed for use with all plate types
- The semi-soft, pliable sealing surfaces prevent the film/foil from being torn or ripped
- Features two sealing surfaces: the shorter length fits inside the walls of raised-rim plates; the longer length is ideal for flat-top and deep-well plates
- Two unique sized corner radii: 1/8" and 1/16" to seal around the plate perimeter; the tighter radius is for use with raised-rim PCR plates
- The non-flat contoured side is designed to press against the plate; the center raised textured surface provides a non-slip finger hold

Bundling Opportunities

- Adhesive films and foils
- Microplates
- PCR plates
- Deep well plates
- Pipettors and pipet tips