

AlumaSeal[®] 96[™]

Easy-to-pierce for direct sample recovery


i Application:

PCR, long-term storage, light sensitive assays and robotics

🌡️ Temperature Range:

-80°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
F-96-100	AlumaSeal 96	100 Sheets



Features and Benefits

- Certified DNase- & RNase- free
- Provides a secure seal to minimize evaporation, reduce cross-contamination and prevent spillage
- Ideal for use with 96-well raised-rim and non-skirted PCR plates
- Features one partial end-tab with no perforation
- Pierceable with a pipet tip or robotic probe for direct sample recovery
- A complimentary sealing paddle is included in each box

Bundling Opportunities

- 96-well PCR plates
- PCR reagents and buffers
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs