

AlumaSeal® CS™


Easy-to-pierce for direct sample recovery

i Application:

Long-term storage & light sensitive assays

🌡️ Temperature Range:

-80°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
FC-100	AlumaSeal CS	100 Sheets
FCS-25	AlumaSeal CS, Sterile	50 Sheets

Features and Benefits

- Certified DNase- & RNase- free
- Specially formulated for cold storage to -80°C
- Excellent barrier properties to retard evaporation
- Easily pierceable with pipet tips or robotic probes for sample recovery
- Not recommended for PCR or thermal cycling

Bundling Opportunities

- Deep well/ storage plates
- Pipettors and pipet tips
- Reservoirs