

# EZcap™ 1x8 PCR FilmStrips™

For sealing a PCR 8-tube strip or a single row of a 96-well plate

**i Application:**  
PCR and short-term storage

**🌡️ Temperature Range:**  
-20°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
US-CS1X8-400	EZcap 1x8 PCR FilmStrips	400 Strips

## Features and Benefits

- Certified DNase- & RNase- free
- Features an extra strong, thick adhesive to minimize evaporation
- Each sheet contains 8 strips on a continuous liner
- An economical and more user-friendly alternative to PCR strip-caps
- Precision cut to a smaller dimension for sealing one PCR 8-tube strip or a single row of a 96-well plate whenever rows must be selectively protected or accessed
- Remove easily with no splash back as can occur with strip-caps

## Bundling Opportunities

- PCR 8-tube strips
- PCR plates
- Deep well/ storage plates
- PCR reagents and buffers
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs