

EZcap™

1x8 PCR

FoilStrips™

For sealing a PCR 8-tube strip or a single row of a 96-well plate

i Application:
PCR and light sensitive assays

🌡️ Temperature Range:
-20°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
F2-CS1X8-300	EZcap 1x8 PCR FoilStrips	300 Strips



Features and Benefits

- Certified DNase- & RNase- free
- Pierceable with a pipet tip or robotic probe for direct sample recovery
- The conformable nature of foil allows for the visual verification of the seal's security around each well
- 2.0mil aluminum foil removes easily without tearing
- Each sheet contains 6 strips on a continuous liner
- An economical and more user-friendly alternative to PCR strip-caps
- Precision cut to a smaller dimension for sealing one PCR 8-tube strip or a single row of a 96-well plate whenever rows must be selectively protected or accessed
- Remove easily with no splash back as can occur with strip-caps

Bundling Opportunities

- PCR 8-tube strips
- PCR plates
- PCR reagents and buffers
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs