

General Purpose Film

For a more secure seal to minimize evaporation

i Application:

General laboratory use, short-term storage, incubation and ELISA

🌡️ Temperature Range:

-40°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
GP-PP-100	General Purpose Film	100 Sheets
GPS-PP-100	General Purpose Film, Sterile	100 Sheets

Features and Benefits

- Certified DNase- & RNase- free
- Features good chemical resistance
- Provides a very secure seal to minimize evaporation, but is slightly harder to remove than SealPlate films
- Less rigid polypropylene film reduces aerosols upon removal
- Specially formulated adhesive minimizes plate residue
- Can be used as a temporary seal while awaiting further processing
- Ideal for use with polypropylene plates

Bundling Opportunities

- ELISA plates
- Plate washers
- Plate readers
- Pipettors and pipet tips
- ELISA reagents and buffers
- Reservoirs