

SealPlate® ColorTab™

For sample identification and plate labeling

i Application:

General laboratory use, short-term storage, incubation and ELISA

🌡️ Temperature Range:

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
SP-IDA-10	SealPlate ColorTab, Assorted	50 Sheets
SP-ID"X"-100*	SealPlate ColorTab, Color	100 Sheets

**Replace the "X" in Catalog # with the specific color ID: (L) Lavender, (B) Blue, (G) Green, (O) Orange, (R) Red*



Features and Benefits

- Provides a secure seal and is easy to remove
- Rigid polyester film features good clarity and easy handling
- Each end-tab features a colored dot and white writing area
- An inscribed color-coded end-tab can be removed and pasted into a laboratory notebook while the other remains on the plate for identification
- Finely perforated end-tabs allow the film sheet to be removed from the plate without disturbing the identification label
- Available in 5 different individual colors or the assortment pack

Bundling Opportunities

- ELISA plates
- Deep well/ storage plates
- ELISA reagents and buffers
- Plate washers
- Plate readers
- Pipettors and pipet tips
- Reservoirs