

SealPlate® RollMate™

Pre-scored film in a convenient, easy-to-use roll format

i Application:

General laboratory use, short-term storage, incubation and ELISA

🌡️ Temperature Range:

-40°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
SP2-RL-500	SealPlate RollMate	500/Roll

Features and Benefits

- Provides a secure seal and is easy to remove
- Rigid polyester film features good clarity and easy handling
- Prescored film on a roll is a more economical and efficient alternative to rapidly seal microplates or racked test tubes
- Faster than applying single sheeted film
- Self-standing, non-slip supports hold the pre-scored film roll securely in place for quick and easy bench-top dispensing
- The roll unwind tension can be easily adjusted to user preference
- Continuous backing liner can be torn off or fed directly into a waste basket for disposal
- Can be used as a temporary seal while awaiting further processing

Bundling Opportunities

- ELISA plates
- ELISA reagents and buffers
- Plate washers
- Plate readers
- Pipettors and pipet tips
- Reservoirs