# ThermalSeal®

An easier to remove film



#### **Application:**

PCR and general storage



### **Temperature Range:**

-40°C to +120°C

CATALOG#	DESCRIPTION	QUANTITY
100-THER-PLT	ThermalSeal	100 Sheets
STR-THER-PLT	ThermalSeal, Sterile	100 Sheets



#### **Features and Benefits**

- · Certified DNase- & RNase- free
- Features a moderate strength adhesive for ease of removal
- A less rigid polypropylene film that minimizes aerosols upon removal
- Features good chemical resistance
- Only recommended for structurally rigid polypropylene PCR plates that don't warp

## **Bundling Opportunities**

- PCR plates
- PCR reagents
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs

