

# eXTReme™ FoilSeal™ Large

For high performance under extreme conditions



## Application:

Long-term storage and light sensitive assays



## Temperature Range:

-80°C to +120°C

| CATALOG #      | DESCRIPTION             | QUANTITY   |
|----------------|-------------------------|------------|
| XTR-FOILLG-100 | eXTReme FoilSeal, Large | 100 Sheets |



## Features and Benefits

- Certified DNase- & RNase- free
- Excellent chemical resistance (e.g. DMSO)
- Recommended for -80°C long-term storage
- Aluminum foil seal offering superior performance in ultra-low temperatures
- 1.9mil thick aluminum foil removes easily without tearing and allows for visual verification of the seal's security around each well
- Specially formulated adhesive provides an excellent seal to minimize evaporation and well-to-well contamination
- Larger cut is designed for greater plate coverage and easier positioning of the film to the plate
- Not recommended for PCR

## Bundling Opportunities

- Storage plates
- Deep-well plates
- Pipettors and pipet tips
- Reservoirs