

eXTReme™ Seal

Extra strong adhesive minimizes evaporation without residue

i Application:
PCR and robotics

🌡 Temperature Range:
-20°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
XTR-100	eXTReme Seal	100 Sheets

Features and Benefits

- Certified DNase- & RNase- free
- Extra-strong thick adhesive minimizes evaporation, yet removes easily without leaving residue on the plate
- Chamfered corners eliminate plate overhang for robotics applications and tight storage conditions
- Ideal for use with raised-rim and non-skirted PCR plates

Bundling Opportunities

- PCR plates
- PCR reagents
- Thermal cyclers
- Pipettors and pipet tips