

ThermalSeal® RT2™

Effective: February 25, 2014
Document No.: MFG-PS-402

Dimensional Information and Tolerances

CRITERIA	DIMENSIONAL INFORMATION / SPECIFICATIONS
Sheeted Film:	
Overall (L x W)	5.625" +/- 0.015 x 3.125" +/- 0.015
End-Tabs	0.375" +/- 0.015
Perforations	0.413" +/- 0.015 from edge of film
Corner Radii	3/32"

Construction

LAYER	DESCRIPTION	THICKNESS
Top Coat	N/A	N/A
Film	Polyester, Clear	2.0 mil +/-10%
Adhesive	Acrylic	1.0 - 1.3 mil
Release Liner	Polyester, Clear	2.0 mil +/-10%

Applications

- qPCR (Real-Time)

Features & Benefits

- "Brilliant" ultra-high optical clarity
- Certified free of detectable levels of DNase-, RNase- and Human/Mouse DNA contamination (PCR Certified)
- Recommended for temperatures from -40°C to 120°C