

# Multi-Swipe™ Sealing Paddle

Effective: February 15, 2022  
Document No.: MFG-PS-433

## Dimensional Information and Tolerances

CRITERIA	DIMENSIONAL INFORMATION / SPECIFICATIONS
Overall Size:	
Overall (L x W)	3.50" ±0.015 x 2.50" ±0.015
Corner Radii	2x 1/16", 2x 1/8"

## Construction

MATERIAL	
	Polyethylene

## Applications

- Microplate sealing

## Features & Benefits

- Reusable, recyclable
- Features semi-soft, pliable sealing surfaces that prevent the film/foil from being torn or ripped
- Designed for use with all plate types