

Plate-Sealing Roller

Effective: February 15, 2022
Document No.: MFG-PS-438

Dimensional Information

COMPONENT	DIMENSIONAL INFORMATION / SPECIFICATIONS
Roller	4.00" wide x 1.125" diameter
Handle	~4.00" long; 1.375" wide at base

Construction

COMPONENT	DESCRIPTION
Roller, Outer	Thermoplastic elastomer
Roller, Inner	Wood core
Handle	Plated steel

Applications

- Plate sealing

Features & Benefits

- Ideal for sealing flat-top plates
- Reusable, autoclavable
- Good chemical resistance