

eXTReme™ 2 Seal

Powerful adhesive removes cleanly with minimal evaporation

i Application:
PCR and robotics

🌡 Temperature Range:
-20°C to +120°C

| CATALOG # | DESCRIPTION | QUANTITY |
|-------------|-----------------------|------------|
| XTR2-100 | eXTReme 2 Seal | 100 Sheets |
| XTR2-LG-100 | eXTReme 2 Seal, Large | 100 Sheets |



Features and Benefits

- Certified DNase- & RNase- free
- Next-generation adhesive formula provides better-than-ever sample protection with improved evaporation rates
- Standard cut features chamfered corners to eliminate plate overhang for robotics applications and tight storage conditions
- Standard cut is ideal for use with raised-rim and non-skirted PCR plates
- Large cut is designed for greater plate coverage and easier positioning of the film to the plate
- Large cut is recommended for use with flat-top, semi-skirted and full-skirted PCR plates as well as deep-well plates
- NOT recommended for qPCR

Bundling Opportunities

- PCR plates
- PCR reagents
- Thermal cyclers
- Pipettors and pipet tips