

PRODUCT SPECIFICATION SHEET

TiterSeal™

Effective: April 20, 2022

Document No.: MFG-PS-442

Dimensional Information and Tolerances

CRITERIA	DIMENSIONAL INFORMATION / SPECIFICATIONS
Sheeted Film:	
Overall (L x W)	5.750" ±0.015 x 3.125" ±0.015
End-Tabs	0.375" ±0.015
Perforations	0.413" from edge of film ±0.015
Corner Radii	1/8"

Construction

LAYER	DESCRIPTION	THICKNESS
Top Coat	N/A	N/A
Film	Polyester	2.0 mil ±10%
Adhesive	Acrylic	0.9 to 1.2 mil
Release Liner	Bleached White Poly-Coated Paper	6.9 mil ±10%

Applications

- Elisa
- Incubation
- General Laboratory Use

Features & Benefits

- Eliminates "edge effect" in sensitive ELISA applications
- Minimizes evaporation, prevents spillage and contamination between wells
- Recommended for temperatures from -40°C to +120°C