

# TiterSeal™

Sample protection you can count on

## **i** Application:

General laboratory use, short-term storage, incubation and ELISA

## **🌡️** Temperature Range:

-40°C to +120°C



CATALOG #	DESCRIPTION	QUANTITY
TTS-100	TiterSeal	100 Sheets

## Features and Benefits

- Minimizes evaporation, prevents well-to-well contamination and spillage
- Eliminates the “edge effect” which can occur with plate lids (e.g. greater evaporation in the outside plate wells)
- Provides a secure seal and is easy to remove
- Rigid polyester film features good clarity and easy handling
- Can be used as a temporary economical seal while awaiting further processing
- Features two end-tabs with perforations

## Bundling Opportunities

- ELISA plates
- Plate washers
- Plate readers
- Pipettors and pipet tips
- ELISA reagents and buffers
- Reservoirs