

# AdvanceSeal™ Heat Seal Peelable PCR Foil

Peelable heat seal for use with single sheet heat seal applicators

**i Application:**  
PCR

**🌡 Temperature Range:**  
+4°C to +94°C

CATALOG #	DESCRIPTION	QUANTITY
HS-FPL-100	AdvanceSeal Heat Seal Peelable PCR Foil	100 Sheets



## Features and Benefits:

- Compatible with single sheet heat seal applicators
- Certified free of detectable levels of DNase, RNase, and Human DNA contamination (PCR Certified)
- Provides excellent evaporation protection
- Peelable foil does not tear upon removal
- Compatible with polypropylene and dual component (polycarbonate/polypropylene) PCR plates
- Backing liner prevents sealing surface contamination (remove prior to sealing)

## Recommended Sealing Parameters

PLATE TYPE	TIME	TEMPERATURE
Polypropylene	2 seconds	140°C
Dual Component PCR (Polycarbonate/Polypropylene)	2 seconds	156°C

*Sealing parameters should be optimized based on plate type to obtain a secure seal.*