


Zone-Free™

For easy access without adhesive fouling of pipet tips or probes

i Application:
Short-term storage and automation

Temperature Range:
 -40°C to +90°C (96-well)
-40°C to +40°C (384-well)



CATALOG #	DESCRIPTION	QUANTITY
ZAF-PE-50	Zone-Free 96 Well	50 Sheets
ZAF-384-50	Zone-Free 384 Well	50 Sheets

Features and Benefits:

- A multilayer film with pierceable, adhesive-free zones to help minimize potential contact between samples and adhesive
- For the 96-well, the white polypropylene sublayer has 96 circular areas void of adhesive
- For the 384-well, the white polyester sublayer has 384 square areas void of adhesive
- Perforated end tabs that aid in alignment
- Well aligned for 96 and 384 multichannel pipets and robotics

Bundling Opportunities:

- 96-well plates
- 384-well plates
- Pipettors
- Pipet tips and robotic probes
- Reservoirs

Excel Scientific
18350 George Blvd.
Victorville, CA 92394 USA

Main: +1 760-246-4545
info@excelscientific.com
excelscientific.com

