


ThermalSeal® 2™ and ThermalSeal® A™

Features a soft, extra thick, conformable adhesive

i Application:
PCR

🌡 Temperature Range:
-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
TS2-100	ThermalSeal 2	100 Sheets
TSA-100	ThermalSeal A	100 Sheets



Features and Benefits:

- Certified DNase- and RNase-free
- Features a soft, extra-thick, conformable adhesive that forms a strong seal around raised wells to minimize evaporation and cross-contamination
- Ideal for use with flat-top, semi-skirted, and fill-skirted PCR plates
- Compatible with hard-shell dual component (e.g. polycarbonate/polypropylene) PCR plates
- Large cut allows for easier positioning of the film to the plate via edge of plate alignment vs. centering

Bundling Opportunities:

- PCR Plates
- PCR Reagents
- Thermal Cyclers
- Pipettors and Pipet Tips
- Reservoirs